

AD-A131 844

THE STUDY OF THIN FILMS ON SEMI-INSULATING GALLIUM  
ARSENIDE BY ELLIPSOMETRY(U) AIR FORCE WRIGHT  
AERONAUTICAL LABS WRIGHT-PATTERSON AFB OH


1/1

UNCLASSIFIED

N T MCDEVITT ET AL. JUN 83 AFWAL-TR-83-4034 F/G 14/2

NL

										END			
										DATE			
										FORMED			
										9 83			
										DTIC			



MICROCOPY RESOLUTION TEST CHART  
NATIONAL BUREAU OF STANDARDS-1963-A

ADA131844

AFWAL-TR-83-4034




THE STUDY OF THIN FILMS ON SEMI-INSULATING GALLIUM  
ARSENIDE BY ELLIPSOMETRY

Neil T. McDevitt  
William L. Baun

Mechanics and Surface Interactions Branch  
Nonmetallic Materials Division

June 1983

Final Report for Period January 1982 to December 1982



Approved for public release; distribution unlimited.

DTIC FILE COPY

MATERIALS LABORATORY  
AIR FORCE WRIGHT AERONAUTICAL LABORATORIES  
AERONAUTICAL SYSTEMS DIVISION  
WRIGHT-PATTERSON AFB, OHIO 45433


83 08 26 049

NOTICE


When Government drawings, specifications, or other data are used for any purpose other than in connection with a definitely related Government procurement operation, the United States Government thereby incurs no responsibility nor any obligation whatsoever; and the fact that the government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data, is not to be regarded by implication or otherwise as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture use, or sell any patented invention that may in any way be related thereto.

This report has been reviewed by the Office of Public Affairs (ASD/PA) and is releasable to the National Technical Information Service (NTIS). At NTIS, it will be available to the general public, including foreign nations.

This technical report has been reviewed and is approved for publication.



NEIL T. MCDEVITT  
Materials Research Engineer  
Mechanics & Surface Interactions Branch



STEPHEN W. TSAI, Chief  
Mechanics and Surface Interactions Branch  
Nonmetallic Materials Division

FOR THE COMMANDER



FRANKLIN D. CHERRY, Chief  
Nonmetallic Materials Division

"If your address has changed, if you wish to be removed from our mailing list, or if the addressee is no longer employed by your organization please notify AFWAL/MLBM, W-PAFB, OH 45433 to help us maintain a current mailing list".

Copies of this report should not be returned unless return is required by security considerations, contractual obligations, or notice on a specific document.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER AFWAL-TR-83-4034	2. GOVT ACCESSION NO. AD-A131844	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) THE STUDY OF THIN FILMS ON SEMI-INSULATING GALLIUM ARSENIDE BY ELLIPSOMETRY		5. TYPE OF REPORT & PERIOD COVERED Internal Jan 82 to Dec 82
		6. PERFORMING ORG. REPORT NUMBER
7. AUTHOR(s) Neil T. McDevitt and William L. Baun		8. CONTRACT OR GRANT NUMBER(s) In-House Report
9. PERFORMING ORGANIZATION NAME AND ADDRESS Materials Laboratory Air Force Systems Command Wright-Patterson AFB, Ohio 45433		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS Project 2303 Task 2303Q108 WUD #50
11. CONTROLLING OFFICE NAME AND ADDRESS Air Force Wright Aeronautical Laboratories Materials Laboratory (AFWAL/MLBM) Wright-Patterson AFB, Ohio 45433		12. REPORT DATE June 1983
		13. NUMBER OF PAGES 34
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		15. SECURITY CLASS. (of this report) Unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Ellipsometry    Gallium Arsenide    Dielectric Films    Epitaxial Films Metal Films		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) A computer-grid procedure is discussed where delta and psi values obtained from a thin filmed surface are used to estimate the optical constants of the film-free surface. Thin films can be measured with reasonable accuracy on gallium arsenide; however, the optical constants of these films cannot be obtained.		

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED


SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

## FOREWORD

This technical report was prepared by N. T. McDevitt and W. L. Baun of the Mechanics and Surface Interactions Branch, Nonmetallic Materials Division, Materials Laboratory, Air Force Wright Aeronautical Laboratories. The work was initiated under Project 2303, "Surface Phenomena" and WUD #50, "Surface and Interface Properties," monitored by Dr. T. W. Haas.

This report covers work performed in-house during the period January 1982 to December 1982.

The authors are grateful to Mr. Gary Griffin for his technical assistance with the computer program.



Accession For	
NRLS CRANI	
100 T-2	
Unannounced	
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A	

# TABLE OF CONTENTS

SECTION	PAGE
I INTRODUCTION	1
II EXPERIMENTAL	4
III PROCEDURE	5
IV RESULTS	7
1. Dielectric Films	7
2. Epitaxial Films	12
3. Metal Films	22
V DISCUSSION AND SUMMARY	27
REFERENCES	28

## LIST OF ILLUSTRATIONS

FIGURE		PAGE
1	Reflection of Polarized Light into its Orthogonal Components	2
2	Computed Delta and Psi Relation for Various Complex Refractive Index Values of GaAs at Zero Film Thickness	8
3	Computed Delta and Thickness Relation for Various Refractive Index Values of $\text{Ga}_2\text{O}_3$ on GaAs	10
4	Computer Generated Delta and Psi Grid for GaAs for Film Thicknesses to $50\text{\AA}$	11
5	Computer Generated Delta and Psi Relation for a $\text{Ga}_2\text{O}_3$ Film Having an Imaginary Part	13
6	Relationship of Delta and Psi for Films of Varying Refractive Index on GaAs for Thicknesses up to $2400\text{\AA}$	14
7	Relationship of Delta and Psi for Dielectric Films on GaAs	15
8	Relationship of Delta and Thickness for a Closed Curve of $\text{Ga}_2\text{O}_3$ on GaAs	21
9	Computed Delta and Psi Relation for Epitaxial Film GaAlAs on GaAs for a Thickness of $2500\text{\AA}$	23
10	Relationship of Delta and Thickness Showing a Dampened Curve for GaAlAs	24
11	Relationship of Delta and Psi for Metal Films on GaAs for a Thickness of $400\text{\AA}$	25
12	Computed Delta and Thickness Relationship for Metal Films Ni, Ge, and Au	26



## LIST OF TABLES

TABLE		PAGE
1	Calculated Values for Delta and Psi for Various Thicknesses of $\text{Si}_3\text{N}_4$ on GaAs	16
2	Calculated Values for Delta and Psi for Various Thicknesses of $\text{Ga}_2\text{O}_3$ on GaAs	17
3	Calculated Values for Delta and Psi for Various Thicknesses of $\text{As}_2\text{O}_3$ on GaAs	18
4	Calculated Values for Delta and Psi for Various Thicknesses of $\text{Al}_2\text{O}_3$ on GaAs	19
5	Calculated Values for Delta and Psi for Various Thicknesses of $\text{SiO}_2$ on GaAs	20

## SECTION I

## INTRODUCTION

Optical methods have for a long time been extensively employed in surface studies and one of the more sensitive techniques used in this field is ellipsometry. Ellipsometry is virtually the only method for the direct determination of the optical constants of a large number of materials, and for the detection and quantitative thickness measurement of films deposited on these materials. The mathematical equations used in ellipsometry were formulated at the end of the last century; however, due to the cumbersome trigonometric equations involved in the analyses of these data, the technique, through the use of computers, has only been utilized in the last decade. This particular study was mainly accomplished through the use of McCrackin's (Reference 1) computer program for ellipsometry.

In principle, ellipsometry involves directing a monochromatic beam of linearly polarized light, at oblique incidence, onto a clean, flat reflecting surface and analyzing the state of polarization of the reflected beam. We can be a little more specific by referring to Figure 1. The plane polarized light has been rotated into s and p components, where the s component vibrates perpendicular to the plane of incidence and the p component parallel to it. The interaction of this light beam with a surface is unique and computation of the differing phase and amplitude of the orthogonal components enables the optical constants of a material to be determined. However, the application of electromagnetic theory to the reflection of light from materials containing free electrons requires the use of a complex refractive index. The free electrons cause an absorption of the incident light and the complex portion of the refractive index is justified by the fact that the imaginary part permits an easier solution to the absorption problem. The complex refractive index  $\bar{n}$  is usually written  $\bar{n} = n - ik$ . Both  $n$  and  $k$  are positive numbers with the negative sign an arbitrary choice for the direction of propagation of the electromagnetic wave.




Figure 1. Reflection of Polarized Light into its Orthogonal Components

Interaction of light with this same surface when it is covered with a continuous, transparent (to the wavelength of light used) isotropic film is also unique, and often allows for the determination of the thickness of the film or its refractive index. In the case of a non-absorbing film or substrate,  $k$  will be zero and the film or substrate's optical constant will be designated only by  $n$ .

Our interest in ellipsometry is aimed at the study of thin films on semiconductors, in particular gallium arsenide (GaAs). Ellipsometry appears quite suitable for the study of these films for several reasons: (1) it is nondestructive by nature; (2) it can be utilized at ambient conditions; and (3) substrates can be studied under realistic processing procedures.

Thin films on semiconducting or semi-insulating GaAs are essential in device and circuit fabrication, particularly in FET (Field Effect Transistors) devices. These films help establish the appropriate properties of the GaAs surface for fabrication purposes. The determination of whether the GaAs surface is clean or contains a film, and the thickness of the film, is an important aspect of this technology. The main emphasis of this report will be on the optical characteristics of the  $\langle 100 \rangle$  surface of the GaAs and the effects dielectrics and metal overlayers have on these properties.

## SECTION II

### EXPERIMENTAL

A Rudolph ellipsometer (Model 43702) was used for this study. The experimental details are described in a previous report (Reference 2). A mercury light source was used (546.1nm) and all measurements were performed at an angle of incidence of 70°. Extinction points were obtained from the polarizer and analyzer settings in Zones 1 and 3. All computations were performed on a PRIME 550+ computer. The program is capable of performing nine different ellipsometric computations. Our main use of the program in this study was centered on the computation of delta and psi for the purpose of studying the refractive index and film thickness of semiconductor materials.

All of the semi-insulating GaAs wafers used in this study were obtained commercially. The wafers were cut from boules grown by the liquid encapsulated Czochralski process. The polished wafers are 50mm in diameter, 0.5mm thick, and are oriented on the <100> plane. No dopants were added intentionally.

## SECTION III

## PROCEDURE

It is necessary to have access to a computer to facilitate the computation of the ellipsometric data. Also it must be remembered that the foundation of ellipsometry is buried deeply in theoretical models. These models require the surface of the substrate to be optically smooth and film free to obtain a true refractive index. The film must be optically isotropic, homogeneous, and transparent to the wavelength of the light source. The light source has to be monochromatic. Other problems that may arise, such as precision of measurement or instrumental errors, can be found in the literature (References 3-6). Consequently, data acquisition from real surfaces still leaves the interpretation aspect of the computed data fairly subjective.

In this report we will be dealing primarily with the ellipsometric parameters delta and psi and their dependence on the values of the refractive index of the substrate ( $n_s$ ), the imaginary part ( $k_s$ ), and the thickness ( $d$ ) of the film. All of these are referenced to an angle of incidence of  $70^\circ$  and 546.1nm incident light. As mentioned previously, the complex refractive index is written

$$\bar{n}_s = n_s - ik_s \quad (1)$$

where  $k_s$  is usually referred to as the extinction coefficient. However, the input into the computer program we are using will not take the value of the extinction coefficient ( $k_s$ ), rather it requires the parameter called the absorption coefficient ( $k_s^*$ ).  $k_s^*$  is related to  $k_s$  by the following:

$$k_s^* = \frac{k_s}{n_s} \quad (2)$$

the complex refractive index may then be written as

$$\bar{n}_s = n_s(1 - ik_s^*) \quad (3)$$

Hereafter, this report will always use  $k_s^*$  when referring to the imaginary part of the refractive index.

Ellipsometry is noted for its sensitivity to changes in the surface of materials. The sensitivity of the technique can be estimated from the following equation that defines the penetration depth of the light as

$$d_p = \frac{\lambda}{4\pi n_s k_s^*} \quad (4)$$

the distance of penetration into the material which is measured in a direction normal to the surface and is dependent on the product of  $n_s k_s^*$ . For light of  $5461\text{\AA}$ , and several values obtained from the literature for  $n_s$  and  $k_s^*$  for gallium arsenide, the calculated penetration depth will be between  $800$  and  $1000\text{\AA}$ .

## SECTION IV

## RESULTS

## 1. DIELECTRIC FILMS

Ellipsometric measurements were performed on five semi-insulating gallium arsenide wafers. These wafers were cut from the same boule. Delta and psi values were obtained from five areas on each two inch wafer. These experimental data points are plotted (solid dots) on the graph (Figure 2). The grid shown in Figure 2 was computer generated. It was formed by using a series of  $n$  and  $k^*$  values for gallium arsenide which were obtained from the literature. Each intersection on the grid represents a delta-psi value that would be obtained from a film-free surface ( $d=0$ ). The spread of the data points is small and indicates the surfaces of the wafers are optically homogeneous. Because of the small spread in the data we can see from the graph that a reasonable average value for these points would be  $\bar{n}_s = 3.98(1-i0.14)$ . However, every gallium arsenide wafer, under ambient conditions, will have a film on the surface. The refractive index obtained from the graph will then be an apparent refractive index and not represent a film-free surface. In order to accurately obtain the thickness of a film on a surface the refractive index of the film-free surface must be known with some accuracy. This usually requires the measurement of the optical constants of a film-free surface while in an ultra-high vacuum environment. Other difficulties involved in this type of measurement are the possible damage to the surface while removing the ambient film while under vacuum, and the presence of the windows of the chamber between the light source, sample, and detector.

The following ellipsometric method is proposed as an easy and quick determination, under ambient conditions, to obtain the optical constants of a film-free surface. However, we are not proposing that this method is capable of predicting the absolute value of a film-free surface, but only a method to obtain the refractive index of a reasonably film-free surface of a particular substrate being studied.






Figure 2. Computed Delta and Psi Relation for Various Complex Refractive Index Values of GaAs at Zero Film Thickness

It has been shown empirically that a nonabsorbing film on an absorbing substrate will lower the value of the real part of the refractive index while increasing the value of the imaginary part. Taking this into consideration and knowing that the data in Figure 2 represents a film covered surface, we can obtain the optical constant of a surface with a thinner film by moving to the right on the grid in Figure 2. Arbitrarily, we chose the next intersection on the grid and read  $\bar{n}_s = 4.04(1-i0.1)$ . We can now generate another set of curves in the following manner. Since  $k_s^*$  is very small compared to  $n_s$ ,  $\delta$  and  $\psi$  values will be insensitive to small changes in  $k_s^*$ . Curves may then be obtained by varying the real part of the refractive index,  $4.04 \pm 3\%$ , and keeping the imaginary part (0.1) constant. The film growth will be 0 to  $50\text{\AA}$  with a refractive index,  $n_f = 1.90$ . With a film of this thickness the refractive index can vary by  $\pm 5\%$  and no error will be introduced into the readings (Figure 3). From the above data a second grid can be constructed as shown in Figure 4. Plotting the experimental data on this graph (solid points) show the GaAs wafers to have a film approximately  $15\text{\AA}$  thick with the optical constants for the film-free surface being  $4.04(1-i0.1)$ . Further indications show in Figure 4 that we have a reasonable value for  $n_s$  which can be seen by the position of the open dots. These data points represent  $\delta$ - $\psi$  values for several of the wafers after they have been chemically etched. The points indicate that some of the original surface film has been removed. These same wafers were subjected to a study by X-ray photoelectron spectroscopy and the corrected cross sections for the  $1s$  line of oxygen and the Auger LMM line of gallium indicating only a small amount of oxygen is present on the surface.

By this comparison between the experimental data points and the calculated  $\delta$ - $\psi$  grid diagrams, formed by varying the real part of the refractive index ( $4.04 \pm 3\%$ ) for a  $50\text{\AA}$  film, we have demonstrated a method that should reveal the optical constants of a film-free substrate that are within the experimental error range of the unique value. The rest of the film data presented in this report will be referenced to the substrate optical constant obtained from Figure 4 where  $\bar{n}_s = 4.04(1-i0.1)$ . The assumption that the GaAs surface has a continuous film under ambient




Figure 3. Computed Delta and Thickness Relation for Various Refractive Index Values of  $\text{Ga}_2\text{O}_3$  on GaAs




Figure 4. Computer Generated Delta and Psi Grid for GaAs for Film Thicknesses to 50 Å

conditions seems valid when compared to the work (Reference 7) reported on GaAs at reduced pressures of oxygen where coverage is obtained around three monolayers. The stoichiometry of the film may be  $\text{As}_2\text{O}_3$ ,  $\text{Ga}_2\text{O}_3$  or  $\text{GaAsO}_4$  but when dealing with a thin film these differences are not crucial (Figure 3) to this procedure. We also assumed that  $k_f^* = 0$ , however, the film may be slightly absorbing and the variations of delta and psi for thin films will not be large enough (Figure 5) to make this procedure invalid.

When a value has been established for a particular surface then delta and psi curves may be calculated for thicker films. Figure 6 shows the effect the refractive index of various films will have on the substrate being studied. A film with an index of 2.50 on this substrate will lose its sensitivity with the psi parameter and will require a fit primarily with the delta parameter. These curves were calculated using  $\bar{n}_s = 4.04(1-i0.1)$ , keeping the real part of the film refractive index constant and  $k_f^* = 0$ , while varying the thickness of the film. Figure 7 shows a similar series of curves for films that are commonly used or found on GaAs. The film was calculated for a film of  $1600\text{\AA}$ . Values for these films and several others are presented in Tables 1-5. The delta and psi values will repeat themselves after an increment of one wavelength of the light source used. This type of curve is considered closed. Comparison of delta with film thickness for a system of  $\text{Ga}_2\text{O}_3$  on GaAs is shown in Figure 8. This figure shows the repeat pattern of a closed curve.

## 2. EPITAXIAL FILMS

Ternary and quaternary compounds are becoming increasingly important as CVD films on GaAs. Problems that are usually encountered in this type of film growth are lattice mismatch and dislocations. The use of ellipsometry in connection with vapor-phase film growth has been studied (Reference 8). Most of these films will have some free electron character, thus  $k_f^* \neq 0$ . When films become slightly absorbing the effect of the substrate on the reflected light becomes less pronounced as the film becomes thicker. The delta and psi curves will no longer be closed




Figure 5. Computer Generated Delta and Psi Relation for a  $\text{Ga}_2\text{O}_3$  Film Having an Imaginary Part




Figure 6. Relationship of Delta and Psi for Films of Varying Refractive Index on GaAs for Thicknesses up to 2400 Å




Figure 7. Relationship of Delta and Psi for Dielectric Films on GaAs



TABLE 1  
CALCULATED VALUES FOR DELTA AND PSI FOR VARIOUS THICKNESSES  
OF  $\text{Si}_3\text{N}_4$  ON GaAs

$\text{Si}_3\text{N}_4/\text{GaAs}$   
NF 1.50  
AS 4.04

0.10

THICK	DFL	PSI	RE. OF. PARALLEL	RE. OF. NORMAL
0.0	166.534	12.202	0.18150	-14.516
50.0	148.954	12.948	0.19295	-37.076
100.0	133.381	14.135	0.21036	-51.455
150.0	120.009	15.758	0.23352	-76.775
200.0	108.592	17.538	0.25840	-97.157
250.0	98.766	19.442	0.28367	-122.645
300.0	90.195	21.424	0.30897	-148.226
350.0	82.606	23.483	0.33259	-172.645
400.0	75.776	25.649	0.35428	-196.774
450.0	69.492	27.992	0.37366	-219.645
500.0	63.508	30.621	0.39052	-241.298
550.0	57.465	33.690	0.40472	-261.678
600.0	51.767	37.362	0.41617	-280.726
650.0	46.390	41.748	0.42485	-298.493
700.0	41.712	46.516	0.43073	-314.922
750.0	37.100	50.407	0.43380	-329.959
800.0	32.672	51.331	0.43408	-343.638
850.0	28.370	48.606	0.43156	-355.970
900.0	24.293	43.979	0.42623	-366.970
950.0	20.380	39.264	0.41811	-376.643
1000.0	16.458	35.189	0.40720	-384.983
1050.0	12.710	31.790	0.39355	-391.983
1100.0	9.050	28.919	0.37720	-397.648
1150.0	5.919	26.407	0.35870	-401.983
1200.0	2.719	24.124	0.33704	-404.983
1250.0	0.717	21.983	0.31377	-406.648
1300.0	-1.384	19.940	0.28899	-406.983
1350.0	-3.853	17.984	0.26347	-405.983
1400.0	-6.801	16.140	0.23832	-403.648
1450.0	-10.259	14.506	0.21514	-400.000
1500.0	-14.137	13.178	0.19603	-395.133
1550.0	-18.511	12.312	0.18340	-389.067
1600.0	-23.279	12.044	0.17962	-381.897
1650.0	-28.590	12.421	0.18521	-373.630
1700.0	-34.212	13.373	0.19917	-364.275
1750.0	-40.194	14.757	0.21922	-353.825
1800.0	-46.414	16.422	0.24291	-342.275
1850.0	-52.851	18.258	0.26823	-329.623
1900.0	-59.343	20.194	0.29364	-315.866
1950.0	-65.878	22.205	0.31824	-301.000
2000.0	-72.456	24.299	0.34117	-285.133
2050.0	-79.077	26.523	0.36201	-268.266
2100.0	-85.740	28.960	0.38045	-250.400
2150.0	-92.445	31.738	0.39630	-231.533
2200.0	-99.192	35.022	0.40944	-211.666
2250.0	-105.980	38.971	0.41983	-190.800
2300.0	-112.809	43.574	0.42744	-168.933
2350.0	-119.679	48.223	0.43224	-146.066
2400.0	-126.589	51.212	0.43424	-122.200

TABLE 2

CALCULATED VALUES FOR DELTA AND PSI FOR VARIOUS THICKNESSES  
OF  $\text{Ga}_2\text{O}_3$  ON GaAs

GA2O3//GAAS

NF 1.00

NS 4.04

0.10

THICK	DEL	PSI	REF. OF PARALLEL	REF. OF NORMAL
0.0	166.534	10.202	0.10148	-14.516
50.0	140.007	12.055	0.19307	-33.534
100.0	133.508	14.202	0.21127	-51.750
150.0	120.227	15.795	0.23406	-66.550
200.0	108.913	17.576	0.25421	-79.901
250.0	99.197	19.487	0.27504	-91.728
300.0	90.742	21.475	0.31033	-102.317
350.0	83.278	23.532	0.33427	-112.579
400.0	76.588	25.689	0.35630	-121.151
450.0	70.470	28.010	0.37608	-129.576
500.0	64.699	30.603	0.39332	-137.617
550.0	59.555	33.617	0.40745	-145.342
600.0	52.711	37.233	0.41984	-152.799
650.0	45.041	41.587	0.42904	-160.002
700.0	34.396	46.506	0.43545	-167.157
750.0	18.844	50.968	0.43911	-174.206
800.0	358.443	52.824	0.44000	178.774
850.0	338.181	50.747	0.43814	171.766
900.0	322.861	46.170	0.43352	164.732
950.0	312.359	41.203	0.42614	157.515
1000.0	304.729	36.829	0.41600	150.211
1050.0	298.452	33.194	0.40313	142.671
1100.0	292.622	30.155	0.38757	134.848
1150.0	286.720	27.532	0.36941	126.661
1200.0	280.427	25.176	0.34882	118.058
1250.0	273.508	22.986	0.32668	108.760
1300.0	265.749	20.902	0.30155	98.747
1350.0	256.912	18.899	0.27599	87.740
1400.0	246.700	16.994	0.25024	75.441
1450.0	234.757	15.242	0.22570	61.484
1500.0	220.728	13.743	0.20424	45.498
1550.0	204.471	12.440	0.18825	27.325
1600.0	186.424	12.084	0.18017	7.398
1650.0	167.813	12.172	0.18153	-13.105
1700.0	150.171	12.885	0.19204	-32.637
1750.0	134.516	14.102	0.20982	-50.207
1800.0	121.084	15.667	0.23237	-65.592
1850.0	109.647	17.447	0.25742	-79.045
1900.0	99.831	19.351	0.28324	-90.546
1950.0	91.297	21.334	0.30860	-101.647
2000.0	83.772	23.386	0.33266	-111.425
2050.0	77.034	25.534	0.35484	-120.491
2100.0	70.883	27.841	0.37476	-129.001
2150.0	65.096	30.410	0.39220	-137.077
2200.0	59.364	33.389	0.40701	-144.813
2250.0	53.179	36.957	0.41912	-152.286
2300.0	45.651	41.259	0.42849	-159.560
2350.0	35.278	46.162	0.43500	-166.693
2400.0	20.114	50.713	0.43894	-173.737

TABLE 3

CALCULATED VALUES FOR DELTA AND PSI FOR VARIOUS THICKNESSES  
OF  $\text{As}_2\text{O}_3$  ON GaAs

AS<sub>2</sub>O<sub>3</sub>/GAAS

NF 1.8°

AS 4.04

0.10

THICK	DEL	PSI	REF. OF. PARALLEL	REF. OF. VOF. PL
0.0	166.534	10.202	0.10198	-14.516
50.0	149.371	12.987	0.19352	-33.620
100.0	134.301	14.259	0.21237	-50.662
150.0	121.500	15.876	0.23547	-65.476
200.0	110.494	17.696	0.26206	-78.355
250.0	101.503	19.627	0.28856	-89.156
300.0	93.592	21.618	0.31545	-99.852
350.0	86.700	23.653	0.34073	-109.048
400.0	80.631	25.750	0.36426	-117.638
450.0	75.228	27.952	0.38568	-125.624
500.0	70.346	30.342	0.40476	-133.172
550.0	65.813	33.040	0.42139	-140.371
600.0	61.372	36.216	0.43549	-147.223
650.0	56.574	40.296	0.44705	-153.597
700.0	50.577	44.894	0.45625	-160.574
750.0	41.760	50.626	0.46252	-166.951
800.0	27.338	56.499	0.46646	-173.218
850.0	5.092	60.076	0.46788	-179.545
900.0	340.326	58.697	0.46679	174.121
950.0	322.891	53.476	0.46318	167.702
1000.0	310.933	47.428	0.45705	161.388
1050.0	303.791	42.089	0.44838	154.168
1100.0	298.491	37.739	0.43714	148.199
1150.0	293.870	34.188	0.42375	141.099
1200.0	289.320	31.229	0.40710	134.141
1250.0	284.504	28.656	0.38833	126.644
1300.0	279.211	26.325	0.36722	119.722
1350.0	273.279	24.138	0.34395	110.063
1400.0	266.547	22.039	0.31887	101.121
1450.0	258.821	19.997	0.29250	91.098
1500.0	249.850	18.023	0.26567	79.932
1550.0	239.305	16.159	0.23931	67.279
1600.0	226.808	14.489	0.21527	52.742
1650.0	212.356	13.142	0.19564	36.022
1700.0	195.141	12.078	0.18300	17.132
1750.0	175.928	12.038	0.17951	-3.018
1800.0	158.955	12.464	0.18588	-22.027
1850.0	142.622	13.477	0.20086	-41.212
1900.0	128.540	14.920	0.22206	-57.280
1950.0	116.641	16.639	0.24697	-71.216
2000.0	106.573	18.515	0.27359	-83.387
2050.0	97.968	20.476	0.30047	-94.182
2100.0	90.522	22.486	0.32649	-103.918
2150.0	84.004	24.544	0.35107	-112.839
2200.0	78.239	26.678	0.37372	-121.126
2250.0	73.078	28.949	0.39416	-128.911
2300.0	68.370	31.454	0.41220	-136.300
2350.0	63.916	34.333	0.42775	-143.371
2400.0	59.393	37.782	0.44076	-150.193

TABLE 4

CALCULATED VALUES FOR DELTA AND PSI FOR VARIOUS THICKNESSES  
OF  $Al_2O_3$  ON GaAs

AL2O3//GaAs

AF 1.55

AC 4.04

0.11

THICK	DEL	PSI	REF. COEFF. 1	REF. COEFF. 2	REF. COEFF. 3	REF. COEFF. 4
0.0	166.534	11.202	0.19198	-14.516	0.84153	177.151
50.0	150.901	13.562	0.19379	-32.196	0.84153	177.151
100.0	137.225	14.191	0.21179	-47.947	0.84153	177.151
150.0	125.671	15.741	0.23499	-61.615	0.84153	177.151
200.0	116.002	17.484	0.26170	-73.477	0.84153	177.151
250.0	107.876	19.327	0.28770	-83.296	0.84153	177.151
300.0	100.987	21.212	0.31475	-92.078	0.84153	177.151
350.0	95.099	23.108	0.34110	-101.330	0.84153	177.151
400.0	90.041	25.007	0.36624	-108.963	0.84153	177.151
450.0	85.693	26.917	0.39981	-116.002	0.84153	177.151
500.0	81.974	28.864	0.41153	-122.643	0.84153	177.151
550.0	78.829	30.842	0.43140	-128.899	0.84153	177.151
600.0	76.216	32.866	0.44930	-134.884	0.84153	177.151
650.0	74.102	34.976	0.46516	-140.599	0.84153	177.151
700.0	72.445	37.251	0.47898	-146.057	0.84153	177.151
750.0	71.174	41.568	0.49081	-151.215	0.84153	177.151
800.0	70.152	45.668	0.50066	-156.075	0.84153	177.151
850.0	69.080	50.861	0.50957	-161.654	0.84153	177.151
900.0	67.267	57.489	0.51457	-166.649	0.84153	177.151
950.0	62.821	65.751	0.51869	-171.592	0.84153	177.151
1000.0	48.741	75.155	0.52095	-176.478	0.84153	177.151
1050.0	357.432	80.754	0.52134	176.641	0.84153	177.151
1100.0	307.836	74.709	0.51902	177.755	0.84153	177.151
1150.0	294.362	65.354	0.51663	179.841	0.84153	177.151
1200.0	290.044	57.080	0.51140	183.876	0.84153	177.151
1250.0	288.236	50.438	0.50440	188.878	0.84153	177.151
1300.0	287.119	45.224	0.49541	193.790	0.84153	177.151
1350.0	286.023	41.093	0.48447	198.437	0.84153	177.151
1400.0	284.656	37.736	0.47153	143.018	0.84153	177.151
1450.0	282.880	34.915	0.45659	137.410	0.84153	177.151
1500.0	280.623	32.452	0.43961	131.571	0.84153	177.151
1550.0	277.941	30.221	0.42063	125.456	0.84153	177.151
1600.0	274.493	28.133	0.39966	119.016	0.84153	177.151
1650.0	270.530	26.126	0.37689	112.147	0.84153	177.151
1700.0	265.883	24.158	0.35242	104.785	0.84153	177.151
1750.0	260.452	22.209	0.32655	96.793	0.84153	177.151
1800.0	254.091	20.275	0.29973	89.093	0.84153	177.151
1850.0	246.596	18.376	0.27261	79.189	0.84153	177.151
1900.0	237.691	16.555	0.24611	67.054	0.84153	177.151
1950.0	227.043	14.886	0.22158	54.245	0.84153	177.151
2000.0	214.340	13.482	0.20082	39.436	0.84153	177.151
2050.0	199.527	12.484	0.18603	22.552	0.84153	177.151
2100.0	183.142	12.030	0.17938	4.120	0.84153	177.151
2150.0	166.430	12.205	0.18202	-14.633	0.84153	177.151
2200.0	150.808	12.969	0.19349	-32.323	0.84153	177.151
2250.0	137.145	14.200	0.21190	-48.036	0.84153	177.151
2300.0	125.604	15.752	0.23506	-61.607	0.84153	177.151
2350.0	115.946	17.495	0.26089	-73.528	0.84153	177.151
2400.0	107.829	19.339	0.28768	-83.858	0.84153	177.151

TABLE 5  
CALCULATED VALUES FOR DELTA AND PSI FOR VARIOUS THICKNESSES  
OF SiO<sub>2</sub> ON GaAs

SiO <sub>2</sub> //GAAS						
ΔE	1.45					
ΔC	4.04					
THICK	DEL	PSI	REF.CF. PARALLEL	REF.CF. PERP.	REF.CF. 45°	REF.CF. 135°
0.0	166.534	12.202	0.18158	-14.716	0.84153	178.950
50.0	152.721	12.201	0.19250	-30.808	0.84051	170.451
100.0	139.796	14.016	0.20919	-45.471	0.83845	174.733
150.0	129.133	15.425	0.23047	-58.293	0.83539	170.483
200.0	120.168	17.519	0.25437	-69.444	0.83100	170.488
250.0	112.625	18.716	0.27905	-79.243	0.82543	168.132
300.0	106.241	20.458	0.30533	-87.058	0.81847	165.803
350.0	100.804	22.211	0.33372	-95.015	0.80996	163.476
400.0	96.158	23.858	0.35533	-103.013	0.79967	161.139
450.0	92.190	25.493	0.37880	-109.642	0.78734	158.168
500.0	88.923	27.424	0.40091	-115.818	0.77164	155.438
550.0	86.008	29.168	0.42150	-121.470	0.75518	152.403
600.0	83.715	30.852	0.44046	-127.197	0.73447	149.080
650.0	81.932	32.817	0.45782	-132.467	0.70994	145.401
700.0	80.662	34.814	0.47349	-137.501	0.68190	141.417
750.0	79.919	37.014	0.48747	-142.351	0.64855	137.489
800.0	79.730	39.513	0.49981	-147.105	0.61003	133.164
850.0	80.135	42.436	0.51053	-151.688	0.55140	128.174
900.0	81.189	45.946	0.51968	-156.141	0.50278	122.850
950.0	82.983	50.252	0.52722	-160.543	0.43847	116.474
1000.0	85.478	55.599	0.53328	-164.853	0.36515	109.468
1050.0	89.655	62.228	0.53782	-169.117	0.28322	101.738
1100.0	96.090	70.549	0.54088	-173.323	0.19421	92.587
1150.0	110.756	79.310	0.54248	-177.514	0.10241	71.731
1200.0	187.410	85.425	0.54262	178.354	0.04342	-9.106
1250.0	249.452	78.088	0.54130	174.116	0.11419	-75.336
1300.0	262.122	69.053	0.53851	169.406	0.20614	-92.016
1350.0	268.070	61.145	0.53425	165.600	0.29438	-102.410
1400.0	271.803	54.625	0.52849	161.362	0.37524	-110.441
1450.0	274.318	49.359	0.52121	156.955	0.44738	-117.523
1500.0	275.947	45.103	0.51238	152.541	0.51053	-123.406
1550.0	276.838	41.616	0.50156	147.980	0.56506	-128.858
1600.0	277.073	38.692	0.48952	143.253	0.61170	-133.780
1650.0	276.704	36.172	0.47625	138.454	0.65137	-138.250
1700.0	275.767	33.936	0.46090	133.437	0.68497	-142.325
1750.0	274.283	31.891	0.44389	128.210	0.71339	-146.772
1800.0	272.260	29.970	0.42521	122.775	0.73738	-149.425
1850.0	269.694	28.122	0.40492	116.946	0.75764	-152.726
1900.0	266.562	26.312	0.38309	110.844	0.77471	-155.718
1950.0	262.821	24.515	0.35986	104.295	0.78907	-158.526
2000.0	258.401	22.720	0.33544	97.224	0.80112	-161.178
2050.0	253.202	20.926	0.31017	89.503	0.81117	-163.699
2100.0	247.075	19.145	0.28449	82.965	0.81947	-166.111
2150.0	239.821	17.409	0.25907	71.390	0.82623	-168.431
2200.0	231.188	15.768	0.23483	60.510	0.83163	-170.678
2250.0	220.895	14.099	0.21302	48.029	0.83578	-172.867
2300.0	208.745	13.105	0.19526	33.723	0.83878	-175.012
2350.0	194.827	12.307	0.18341	17.780	0.84070	-177.126
2400.0	179.758	12.012	0.17907	0.534	0.84159	-179.224




Figure 8. Relationship of Delta and Thickness for a Closed Curve of  $\text{Ga}_2\text{O}_3$  on GaAs

for this type of film. The delta and psi values will approach the values representative of the bulk film with increasing film thickness. This type of curve, shown in Figure 9, is calculated for GaAlAs/GaAs where  $\bar{n}_f = 4.2(1-i0.067)$ . A comparison of delta and thickness values for this system (Figure 10) shows a dampened curve as the substrate becomes increasingly obscured by the film. The magnitude of the absorption character of the film would determine the usefulness of ellipsometry for each epitaxial film studied. For the film data shown in Figure 10 information could be reasonably obtained from a 200Å film and possibly up to 600Å. Thicker films would become increasingly difficult to interpret using null ellipsometry.

### 3. METAL FILMS

The use of ellipsometry in connection with programs dealing with contacts and interconnects on compound semiconductors was reviewed. A short study was performed looking at a small number of metal films on GaAs. The absorption character ( $k_f^*$ ) of metallic films is usually greater than 0.25. The limitation of the sensitivity of ellipsometry will depend primarily on the value of  $k_f^*$  of each film studied. Knowing this limitation will be an important aspect of a study concerning metal films. The metal films looked at in this study were, nickel 1.4-(1-i1.8), gold 0.43(1-i5.12), and germanium 5.46(1-i0.32). The delta and psi curves (Figure 11) are characteristic of metal films and show an entirely different response when compared to dielectric and epitaxial films. These data were generated for 400Å thick films and a comparison of delta and thickness is reported in Figure 12. Sensitivity for the GaAs surface will become obscure for nickel and gold around 150Å, while germanium films may be studied to 300Å. However, the germanium film will suffer in sensitivity where the maximum occurs in the curve.




Figure 9. Computed Delta and Psi Relation for Epitaxial Film GaAlAs on GaAs for a Thickness of 2500Å






Figure 10. Relationship of Delta and Thickness Showing a Dampened Curve for GaAlAs




Figure 11. Relationship of Delta and Psi for Metal Films on GaAs for a Thickness of 400Å




Figure 12. Computed Delta and Thickness Relationship for Metal Films Ni, Ge, and Au

## SECTION V

### DISCUSSION AND SUMMARY

Ellipsometric data has been obtained from a number of commercially prepared gallium arsenide wafers. The wafers were 50mm in diameter and had a polished surface. The optical constants were measured from the as-received surfaces. A grid procedure consisting of the as-received delta and psi readings and literature values for  $n_s$  and  $k_s^*$  is proposed for finding the optical constants of a film-free surface.

Using the optical constant for the film-free surface a series of experimental and computer studies were performed for dielectric, epitaxial, and metal films on gallium arsenide. When dealing with very thin dielectric films, good fits can be obtained between observed and calculated data giving a reasonable measure of film thickness. However, the optical constants for these thin films cannot be obtained with any accuracy by ellipsometry, because delta and psi approach the same values regardless of the optical constants of the film as the thickness tends towards zero.

# REFERENCES

1. F. McCrackin, A Fortran Program for Analysis of Ellipsometer Measurements, Natl. Bur. of Std. Technical Note 479, U.S. Government Printing Office, Washington D. C. (1969).
2. N. T. McDevitt, Air Force Materials Laboratory, Technical Report AFML-TR-73-245, January 1974, Wright-Patterson AFB, Ohio.
3. E. Schmidt, J. Opt. Soc. Am. 60, 490(1970).
4. F. L. McCrackin, J. Opt. Soc, Am. 60, 57(1970).
5. R. J. Archer and C. V. Shank, J. Opt. Soc, Am. 57, 191(1967).
6. H. T. Yolken, R. M. Waxler, and J. Kryger, J. Opt. Soc. Am 57, 283(1967).
7. I. Lindau and W. E. Spicer, in "Electron Spectroscopy: Theory, Techniques and Applications," C. Brundle and A. Baker, Eds., Vol. 4, Chapter 4, Academic Press, New York 1981.
8. J. B. Theeten, Surface Sci. 96, 275(1980).

ATE  
LME